

## CT-6 Leg Splint Visual User Guide



FareTec Inc  
1610 West Jackson St  
Painesville OH 44077  
440 350 9510 Voice  
440 350 9520 Fax  
800 322 3273 Toll Free

Copyright ©2010 Fare Tec. Inc.

### Assess Injury

#### Step 1


Check circulation



Cut off pant leg if needed.



Assess injury to determine if traction is needed.  
Determine pain level. Compare lengths of the patient's legs.




### Assemble

#### Step 2

Remove splint from bag.



### Components



#### Repacking:

1. Slide leg straps to end of splint, but not over tube edges.
2. Pull apart tube sections and fold each on top of the other. Allow straps to hang.
3. Fold ankle hitch and ischial strap alongside folded tube sections.
4. Roll up splint tightly with leg straps and slide into bag, ankle hitch first.

#### Care and Use:

USE PPE : Gloves, mask and eye protection.  
CLEAN : Wash with soap and water first and dry with disposable towels.  
DISINFECT : Use commercial disinfectant according to directions or solution of 1 part bleach to 100 parts water. Treat entire CT-6 and dry with disposable towels.  
DISCARD : Dispose of used cleaning materials and all PPE (except eye protection) in properly

### Assemble

#### Step 3


Shake splint up and down to connect the tube sections.  
Manually connect any unlinked sections.



### Sizing

#### Step 4

Place unit alongside uninjured leg. Splint should extend just above hip crest and 6" in. below patient's foot.



Two straps above knee. Two straps below knee.