12 V-Stecker / 12 V plug / Prise 12 V

für / for / pour

Ladeschnittstelle BASE STATION

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BASE STATION charger interface

ACCUVAC connecting cable

Interface de charge BASE STATION

Câble de raccordement ACCUVAC

Montageanleitung
Assembly instructions
Notice de montage

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These assembly instructions describe how to fit an alternative 12 V plug to the following charging interface/connecting cable:

WM number	Charging interface/connecting cable for
WM 22436	BASE-STATION 1 NG BASE-STATION 3 NG BASE-STATION II BASE-STATION mini II
WM 10650	ACCUVAC Rescue ACCUVAC Basic

Prerequisites for an alternative 12 V plug

- Suitable for a 13.8 V on-board power supply
- The temperature range must correspond to the internal temperatures arising during use
- Current load > 6 A
- Strain relief is essential

Remove existing 12 V plug

1. Cut connecting cable approx. 1 cm behind the 12 V plug.
2. Dispose of the plug you have cut off as electronic waste.


Attach the alternative 12 V plug

1. Remove the insulation from the cable ends.
2. Twist cable shielding and black cable together, then crimp with a cable end-sleeve.
3. Twist red and brown cable together, then crimp with a cable end-sleeve.

4. Connect the black cable/cable shielding to the negative contact and the red-and-brown cable to the positive contact of the alternative 12 V plug.

Test the alternative 12 V plug

1. Measure the electrical passage of the following connections:

- Positive contact of connecting block/connector against positive contact of alternative 12 V plug
 - Negative contact of connecting block/connector against negative contact of alternative 12 V plug
2. Measure both contacts on the BASE STATION connecting block/ ACCUVAC connector for short-circuiting. There may not be any short circuit.
 3. Connect the alternative 12 V plug to the vehicle's on-board power supply.
 4. Check the polarity of the voltage at the contacts of the BASE STATION connecting block/ACCUVAC connector.

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GmbH+Co. KG**